



Chiropractic Radiographic System with Anatomical Programming

High-frequency X-ray system built for the chiropractic practice

Model: DC-1



A chiropractic X-ray room with a high-frequency generator rated at 500 mA and 600 mAs. High-frequency output lowers patient dose, sharpens images and helps the tube last longer. Technique is set by anatomy in three steps, and Sharp Imaging Solutions can convert the room to Canon CXDI DR at any time.

■ Three-step setup

Select the anatomical region, select the view, enter the patient measurement. The display shows the optimum mA, kVp and mAs.

■ Reliable

Built-in self-checks monitor each exposure as it happens.

■ Built to last

Tube stands and vertical structures are solid steel, with stainless steel bearings for accurate alignment of the X-ray tube and bucky.

■ Ready for DR

Compatible with most digital imaging systems. Add a Canon CXDI DR upgrade from Sharp Imaging Solutions at any time.

■ Warranty

7-year parts warranty from the manufacturer on film-based systems.

Specifications



Generator console with anatomical programming

System	
Generator	HF-30B, high frequency, 500 mA, 600 mAs output
Technique selection	Anatomical programming with researched and proven techniques
Structure	Solid steel tube stand and vertical cassette holder, stainless steel bearings
Options	Traditional wall stand
Digital	Upgradeable to DICOM-compliant DR

Sold, installed and serviced by Sharp Imaging Solutions

- **Site planning and installation.** Room review, deinstallation of old equipment, installation and calibration.
- **Canon CXDI DR.** Authorized Canon dealer and service provider for detectors, software and DR upgrades.
- **Connectivity.** DICOM worklist and PACS setup, with training for your technologists.
- **Service and parts.** Preventive maintenance, repairs and parts for multi-vendor X-ray equipment.

Request a quote

561-441-2157

sales@sharpimagingsolutions.com
sharpimagingsolutions.com

Manufactured by Bennett X-Ray Systems. Sharp Imaging Solutions LLC is an independent dealer that sells, installs and services this equipment, and is not the manufacturer.

Specifications are based on the manufacturer's published information and may change without notice. Final configuration, specifications and warranty terms are confirmed in your written quote.

Bennett and DC-1 are trademarks of their owner. Canon, CXDI and related names are trademarks of Canon Inc. Windows is a trademark of Microsoft Corporation. These names are used only to identify the equipment and do not imply endorsement.

Product photos are the property of the equipment manufacturer and are used with permission.

© 2026 Sharp Imaging Solutions LLC. Brochure text and layout by Sharp Imaging Solutions LLC.