



Floor-Mounted Tube Stand Radiographic Room

Full capability in rooms with low ceilings

Model: **FMT Classic**



A cost-effective room for busy clinics with limited space. It needs a ceiling height of just 92.5 in., so most existing rooms work without construction, and Sharp Imaging Solutions can add Canon CXDI DR now or later.

■ Table, wall stand or wheelchair

The tube stand images on the table or wall stand and accommodates wheelchair exams.

■ Quick positioning

40-in. and 72-in. electronic detents, a backlit angulation dial and a one-touch release handle.

■ Patient-friendly

Elevating four-way float table rated for 400 lb, recessed foot pedals and a comfortable chin rest.

■ Low cost of ownership

Backed by a 5-year parts warranty from the manufacturer.



Tube stand controls and angulation dial



High-frequency generator and console

Specifications

Radiographic table (S211)		Floor-mounted tube stand (J700)	
Tabletop	84 x 30 in, flat and seamless, four-way float	Rails	10 ft
Height	23 in. to 35 in.	Rotation	360° base with 90° detents
Travel	±5 in. transverse, 31 in. longitudinal	Transverse tube arm	8.5 in. travel
Controls	Foot treadles for tabletop and height	Controls	Angulation dial, operator handgrips, electric locks, platform tube mount
Bucky	Grid cabinet with electric lock; 10:1, 103 LPI grid	Wall stand (S109)	
Patient load	400 lb	Type	Wall-mounted, fully counterbalanced, center mount
Power	120 VAC	Bucky travel	17 x 17 in. grid cabinet, floor to 76 in.
		Grid	10:1, 103 LPI, specify focal range
		Locks	Electric locks, left or right load, 24 VDC
		Generator	
		• 40 kW to 80 kW high-frequency generator options	

Sold, installed and serviced by Sharp Imaging Solutions

- **Site planning and installation.** Room review, deinstallation of old equipment, installation and calibration.
- **Canon CXDI DR.** Authorized Canon dealer and service provider for detectors, software and DR upgrades.
- **Connectivity.** DICOM worklist and PACS setup, with training for your technologists.
- **Service and parts.** Preventive maintenance, repairs and parts for multi-vendor X-ray equipment.

Request a quote

561-441-2157

sales@sharpimagingolutions.com
sharpimagingolutions.com

Manufactured by AMRAD Medical, a brand of Summit Industries, LLC. Sharp Imaging Solutions LLC is an independent dealer that sells, installs and services this equipment, and is not the manufacturer.

Specifications are based on the manufacturer's published information and may change without notice. Final configuration, specifications and warranty terms are confirmed in your written quote.

AMRAD Medical, FMT Classic and related product names are trademarks of Summit Industries, LLC. Canon, CXDI and related names are trademarks of Canon Inc. Windows is a trademark of Microsoft Corporation. These names are used only to identify the equipment and do not imply endorsement.

Product photos are the property of the equipment manufacturer and are used with permission.

© 2026 Sharp Imaging Solutions LLC. Brochure text and layout by Sharp Imaging Solutions LLC.