



Mobile C-Arm, All-in-One

Compact flat-panel fluoroscopy for tight spaces

Model: **Agility**



A single-cabinet flat-panel C-arm that moves easily between rooms. Images appear on a 27-in. 4 MP monitor mounted on an adjustable arm, so the surgical team can see them from anywhere around the table.

■ Flat panel detector

Acquisition areas from 21 x 21 cm to 30 x 30 cm, with a removable anti-scatter grid for low-dose work.

■ Pulsed, modulated power

Pulsed output adjusts to the procedure, keeping dose low without giving up image quality.

■ Easy access


The touch console swivels, so it can be used from any side of the C-arm while the live image stays in view.

■ Remote display

Wi-Fi video to an optional wireless monitor cart, ceiling suspension or wall monitors.



Optional wireless monitor cart



In the operating room

Specifications

Movement

Orbital rotation	160° (+65° to -95°)
Angular rotation	More than ±210°
Vertical travel	18 in. (450 mm)

Display and console

Main monitor	27 in., 4 MP, live and stored images side by side; height, angle and tilt adjustable
Console	High-brightness touchscreen for selecting, editing, sending and zooming images
Wireless monitor cart	Optional; lithium-ion battery, duplicates the main monitor

Quality features

- Laser localizers on the tube and detector for zero-dose centering
- Removable grid for low-dose pediatric and gridless procedures
- Foot switch or multi-function hand switch; optional wireless footswitch
- USB, Ethernet, HDMI and Wi-Fi connectivity; optional integrated printer
- Medical-grade certified display; instant login with optional NFC
- Native remote connectivity for efficiency tests, software updates and maintenance

Sold, installed and serviced by Sharp Imaging Solutions

- **Site planning and installation.** Room review, deinstallation of old equipment, installation and calibration.
- **Canon CXDI DR.** Authorized Canon dealer and service provider for detectors, software and DR upgrades.
- **Connectivity.** DICOM worklist and PACS setup, with training for your technologists.
- **Service and parts.** Preventive maintenance, repairs and parts for multi-vendor X-ray equipment.

Request a quote

561-441-2157

sales@sharpimagingsolutions.com
sharpimagingsolutions.com

Manufactured by AMRAD Medical, a brand of Summit Industries, LLC. Sharp Imaging Solutions LLC is an independent dealer that sells, installs and services this equipment, and is not the manufacturer.

Specifications are based on the manufacturer's published information and may change without notice. Final configuration, specifications and warranty terms are confirmed in your written quote.

AMRAD Medical, Agility and related product names are trademarks of Summit Industries, LLC. Canon, CXDI and related names are trademarks of Canon Inc. Windows is a trademark of Microsoft Corporation. These names are used only to identify the equipment and do not imply endorsement.

Product photos are the property of the equipment manufacturer and are used with permission.

© 2026 Sharp Imaging Solutions LLC. Brochure text and layout by Sharp Imaging Solutions LLC.