



Wall-Mounted Chest Stand System

For pulmonology, occupational medicine and upright exams

Model: CS Elite



A compact room for chest and other upright exams, made up of a tube stand and a counterbalanced wall stand. Electric locks and a collimator light switch at the wall stand keep positioning quick. Add a mobile table if you need table exams too, and add Canon CXDI DR whenever you are ready.

- **Built for upright exams**

Chest and other standing exams, with a choice of 30 kW to 50 kW generators.

- **Space-saving**

Center-mounted or flag-mounted wall stand, floor-to-wall or floor-to-ceiling tube stand.

- **Fast and efficient**

Electric locks on every motion, angulation dial with handgrips, and a remote collimator light.

- **DR upgradeable**

Add a Canon CXDI detector now or later.

Specifications

Tube stand		Wall stand	
Mounting	Floor-to-wall or floor-to-ceiling, non-rotating, with transverse travel	Design	Center-mounted or flag-mounted, fully counterbalanced
Tube arm	24 in., angulating, platform tube mount	Bucky	17 in. grid cabinet with front panel; heavy-duty cassette tray
Transverse travel	±5 in.	Grid	10:1, 103 LPI, specify focal range
Rails	8 ft (6.5 ft of travel)	Features	Electric locks, PA grips, remote collimator light, 24 VDC
Ceiling height	88 in. to 121 in. for floor-to-ceiling mount	Generator 30 kW to 50 kW high-frequency generators	
Vertical travel	13.25 in. minimum (focal spot to floor), 73.5 in. maximum		
Locks	Vertical, horizontal, transverse and angulation electric locks; 24 VDC		

Sold, installed and serviced by Sharp Imaging Solutions

- **Site planning and installation.** Room review, deinstallation of old equipment, installation and calibration.
- **Canon CXDI DR.** Authorized Canon dealer and service provider for detectors, software and DR upgrades.
- **Connectivity.** DICOM worklist and PACS setup, with training for your technologists.
- **Service and parts.** Preventive maintenance, repairs and parts for multi-vendor X-ray equipment.

Request a quote

561-441-2157

sales@sharpimagingolutions.com

sharpimagingolutions.com

Manufactured by AMRAD Medical, a brand of Summit Industries, LLC. Sharp Imaging Solutions LLC is an independent dealer that sells, installs and services this equipment, and is not the manufacturer.

Specifications are based on the manufacturer's published information and may change without notice. Final configuration, specifications and warranty terms are confirmed in your written quote.

AMRAD Medical, CS Elite and related product names are trademarks of Summit Industries, LLC. Canon, CXDI and related names are trademarks of Canon Inc. Windows is a trademark of Microsoft Corporation. These names are used only to identify the equipment and do not imply endorsement.

Product photos are the property of the equipment manufacturer and are used with permission.

© 2026 Sharp Imaging Solutions LLC. Brochure text and layout by Sharp Imaging Solutions LLC.